

 STAARSURGICAL™

EVO ICL™

Patient Selection




# Patient Selection

**The EVO ICL™ family of lenses\* (EVO) is indicated for use in phakic eye treatment in patients 21 – 60 years of age who meet the criteria listed below:**

- The correction/reduction of myopia with or without astigmatism.
- Requiring a lens with sphere power of -0.5 D to -20.0 D and cylinder power up to +6.0 D.
- With an anterior chamber depth (ACD) equal to or greater than 2.8 mm, as measured from the corneal endothelium to the anterior lens capsule.

## How does it work?

EVO is placed directly behind the iris and in front of the natural crystalline lens. In this position, EVO helps the eye to focus light properly onto the retina to create clear distance vision.



## **EVO candidates may include:**

- Myopes starting as low as -0.5 D to -20.0 D.
- Patients being considered for laser vision correction procedures.
- Patients with thinner corneas <sup>1</sup>.
- Patients with dry eye risk factors <sup>2</sup>.
- Patients whose corneal topography is less suited for laser vision correction.

## **Patients not suitable for EVO include those:**

- With an ACD of < 2.8 mm.
- With anterior chamber angle < Grade III as determined by gonioscopic examination.
- Who are pregnant or nursing.
- With Keratoconus.
- Who are less than 21 years of age.
- With low/abnormal corneal endothelial cell density, Fuch's dystrophy or other corneal pathology.
- With any cataract in the operative eye or nontraumatic cataract in the fellow eye.
- With primary open angle or narrow angle glaucoma.

# EVO ICL™

## Patient Selection


- With previous or pre-existing ocular disease that would preclude post operative visual acuity of 0.477 (logMAR), 20/60 (Snellen), 0.33 (Decimal), or better.
- Who are amblyopic or blind in the fellow eye.
- With ocular hypertension in either eye.

\* EVO ICL family of lenses include EVO ICL, EVO+ ICL, EVO Toric ICL, EVO+ Toric ICL

## Benefits



Made with **Collamer™** **biocompatible material** that works in harmony with the natural eye.



**Permanent yet removable** by a doctor, if necessary.




**Excellent vision** day and night.<sup>3,4</sup>

# 99.4 %

of patients surveyed **would have the procedure again.**<sup>5</sup>



**Does not induce** dry eye syndrome.<sup>2</sup>



**UV** protection.



## Preoperative information

### Patient work up

- A standard, full ophthalmic exam should be performed.
- Measurements needed for performing EVO calculations in the Online Calculation and Ordering System (OCOS).

### Measurements recommended For patient assessment and records

- Corneal Endothelial Cell Density (ECD) assessment.
- Gonioscopic assessment of the angle, Grade III or higher.
- Axial length.
- Accurate and stable refraction.

#### Example of OCOS data entry tab

Back vertex distance

Manifest and/or cycloplegic refraction

Keratometry

Anterior chamber depth

Corneal thickness

White to white

Contact lens over refraction sphere (optional)

Calculate for	
<input type="radio"/> ICL	<input checked="" type="radio"/> Toric ICL
Patient ID	123456789
Patient name	John Doe
Operative eye	<input type="radio"/> OD <input checked="" type="radio"/> OS
DOB	1975.01.01
Gender	<input type="radio"/> M <input checked="" type="radio"/> F
BVD	1.2
Sphere	-5.5
Cylinder	2
Axis	90
	Power Degrees
K1	42 @ 0
K2	44 @ 90
ACD	3.0
CT	.5
WW	11.8
CL sphere	0
Any previous intervention?	<input checked="" type="radio"/> No <input type="radio"/> Yes

# Patient Selection

## Postoperative information

### Recommended patient

### Postoperative assessment <sup>6</sup>


Postoperative same day, day 1, day 7 and beyond.

- Visual acuity.
- Intraocular pressure.
- Assess the ICL to crystalline lens vault.
- Biomicroscopy to assess:
  - EVO centration.
  - Inflammation.


### Measuring the vault

Although the postoperative vault of EVO is intended to be approximately equal to the central corneal thickness, the optimal vault should be between 50% and 150% of central corneal thickness, this being equivalent to a range of 250 to 900 microns.


However, in the absence of symptoms, lens vault outside this range may not necessarily require exchange or removal.



Normal vault ◀



Shallow vault ◀



▶ High vault

## EVO models available:

Models	Spherical Power (D)	Cylindrical Power (D) (For EVO/EVO+ Toric)	Overall Diameters (mm)
EVO+	-0.5 to -14.0	0.5 to 6.0	12.1
			12.6
EVO	-0.5 to -18.0		13.2
			13.7


- Spherical models available in 0.25 D increments from -0.5 D to -3.0 D and 0.50 D increments from -3.0 D to -20.0 D
- Toric models are only available in 0.5 D increments

## References

1. Albo C, Nasser T, Szykarski DT, Nguyen N, Mueller B, Libfraind L, Parkhurst G. A Comprehensive Retrospective Analysis of EVO/EVO+ Implantable Collamer Lens: Evaluating Refractive Outcomes in the Largest Single Center Study of ICL Patients in the United States. Clin Ophthalmol. 2024 Jan 9;18:69-78.
2. Ganesh S, Brar S, Pawar A. Matched population comparison of visual outcomes and patient satisfaction between 3 modalities for the correction of low to moderate myopic astigmatism. Clin Ophthalmol. 2017;11:1253-63.
3. Parkhurst GD. A prospective comparison of phakic collamer lenses and wavefront-optimized laser-assisted in situ keratomileusis for correction of myopia. Clin Ophthalmol. 2016;10:1209-15.
4. Martínez-Plaza E, López-Miguel A, López-De La Rosa A, et al. Effect of the EVO+ Visian Phakic Implantable Collamer Lens on Visual Performance and Quality of Vision and Life, Am J Ophthalmol 2021;226: 117–125.
5. Packer M. The Implantable Collamer Lens with a central port: review of the literature. Clin ophthalmology. 2018;12:2427-38.
6. Chuck RS, Jacobs DS, Lee JK, Afshari NA, Vitale S, Shen TT, et al. Refractive Errors & Refractive Surgery Preferred Practice Pattern®. Ophthalmology. 2018;125(1):P1-p104.

### Important Safety Information for EVO/EVO+ ICL

The EVO/EVO+ ICL is indicated for phakic patients 21-60 years of age to correct/reduce myopia up to -20.0 D with up to 6.0 D of astigmatism. Careful preoperative evaluation and sound clinical judgment should be used by the surgeon to decide the risk/ benefit ratio before implanting a lens in a patient with any of the conditions described in the DFU. Prior to surgery, physicians should inform prospective patients of possible risks and benefits associated with the EVO/EVO+ ICL. Reference the EVO/EVO+ ICL DFU available at <https://edfu.staar.com/edfu/> for a complete listing of indications, contraindications, warnings and precautions.



discovericl.com  
©2024 STAAR® Surgical Company  
ES-EVO ICL-25-0002

Hauptstrasse 104  
CH -2560 Nidau, Switzerland  
+41 32 332 88 88



EVO ICL<sup>TM</sup>

 STAAR SURGICAL<sup>TM</sup>